

INSTRUCTIONS

#21 STATIC PROGRESSIVE POSITIONING TO EXTEND MP JOINTS

Provides static progressive positioning to help correct elastic flexion contractures and joint tightness.

Contra indicated for Fixed Contractures, Fragile Skin, Severe Edema, and Unhealed Fractures.

*Tightening top strap will exert a force to extend PIP Joint. Range of motion should be increased very gradually, and tolerated by the patient. Excessive force can cause joint damage, skin breakdown, and pain. Do not be overly aggressive. Depending on condition of patient, establish an appropriate schedule for wearing, skin examination and amount of force.

APPLICATION—Loosen straps. Slide hand into splint. With fingers comfortably resting on platform adjust finger separators. Position and secure strap for fingers to liner on underside of platform. Comfortably wrap strap around fingers and secure to hook on end of strap. Wrap strap around wrist & comfortably secure. Tighten top strap to extend MP joints. Adjust as needed.

ADJUSTMENTS

—“Non slip” pads—

Under fingers—position & attach to liner

On strap for fingers—position & attach to loop on strap

—Extension wedge—


Reflect liner—Attach to platform using hook & loop

Trim as needed—Replace liner

CLEANING

Splint—Clean with mild soap and damp towel. Use damp towel to remove excess soap. Air dry completely before reapplying.

Liner & Pads—Remove from splint. Hand wash cold water and mild soap. Thoroughly rinse. Air dry completely.



◆ Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritation, or have any questions or concerns.

Single patient use. To be used only when ordered by a physician and applied by qualified medical professional

