


INSTRUCTIONS #9 Reverse Knuckle Bender Splint

Open wrist strap and slide hand into splint. Position M.P. extension bar under fingers between M.P. and P.I.P. joints. Position round pad in palm. Bring strap over top of wrist and secure. Finish attaching rubber bands to hooks. Splint is designed to dynamically extend M.P. joints while allowing active flexion. Evaluate fit, function and instruct patient on use.

Modifications: Remove from patient. Rubber bands can be added or subtracted. Wire and band can be carefully bent to improve fit and function of splint.

Clean with mild soap and damp towel. Air dry completely before reapplying.

Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritations or have any questions or concerns.



BUNNELL #9 Reverse Knuckle Bender Splint to extend the metacarpophalangeal joints.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use