

TIBURON ENTERPRISES, INC.
DBA WENIGER
EXCLUSIVE
MANUFACTURER OF
ESTABLISHED 1907

BUNNELL




WRIST-HAND-FINGER ORTHOSIS
STATIC-DYNAMIC – ADJUSTING STATIC POSITIONING
LATEX FREE-SOFT STRAPS-ORTHOSIS FITS RIGHT OR LEFT HAND



TIBURON MEDICAL ENTERPRISES, INC.
DBA WENIGER, MFG OF BUNNELL SPLINTS
Expanded product line— Customer service
Please visit BunnellSplints.com


Sizing information on back page

BUNNELL ACTIVE HAND AND FINGER SPLINTS



#6

Spring Wire Safety Pin Splint to extend PIP Joint. Positional strap and pad.
Sizes: XS, SM, MD, LG, XL




#6A

Modified Safety Pin Splint to extend PIP Joint. Dynamic Spring. Positional strap and pad. Sizes: XS, SM, MD, LG, XL



#6MINI

Mini Safety Pin Splint to extend DIP Joint. Positional strap and pad.
Sizes: 1, 2, 3, 4




#6A MINI

Mini Modified Safety Pin Splint to extend DIP Joint. Dynamic Spring. Positional strap and pad.
Sizes: 1, 2, 3, 4




#5

Finger Extension Splint w/Clockspring to extend a finger or thumb.
Sizes: XS, SM, MD, LG, XL




#5A

Finger Extension Splint w/Wrist Support. Dynamic Spring to M.P. joint for improved function. Malleable support extends to M.P. joint.
Sizes: XS, SM, MD, LG, XL




#11
Finger Knuckle Bender Splint to flex the proximal interphalangeal joint.
Sizes: XS, SM, MD, LG, XL



#11A
Reverse Finger Knuckle Bender Splint to extend the proximal interphalangeal joint.
Sizes: XS, SM, MD, LG, XL



#2
Knuckle Bender Splint to flex the metacarpophalangeal joints.
Sizes: XS, SM, MD, LG, XL




#2A
Knuckle Bender Splint to flex the metacarpophalangeal joints w/outrigger to extend the distal joints.
Sizes: XS, SM, MD, LG, XL




#2B
Knuckle Bender Splint to flex the metacarpophalangeal joints with extra segment to flex the distal joints.
Sizes: XS, SM, MD, LG, XL



#10
Composite Spring and Elastic Splint to extend the distal interphalangeal joints.
Sizes: SM, MD, LG




#9
Reverse Knuckle Bender Splint to extend the metacarpophalangeal joints.
Sizes: XS, MD, LG, XL



#9A
Reverse Knuckle Bender Splint to extend the metacarpophalangeal joints with outrigger to extend the distal joints.
Sizes: XS, SM, MD, LG, XL




#1
Spring Cock-up
Splint to dorsiflex the wrist.
Kydex plastic forearm shell.
Sizes: SM, MD, LG




#1A
Spring Cock-up
Splint to dorsiflex the wrist with
outrigger to flex the M.P. joints.
Kydex plastic forearm shell.
Sizes: SM, MD, LG




#8
Dorsal Wrist
Splint to hold the wrist in flexion
after repair of flexor tendons or
after median and ulnar nerve
suture. Kydex plastic forearm
and wrist section.
Sizes: SM, MD, LG



#8A
Dorsal Wrist
Splint to hold the wrist in flexion
with outrigger to extend distal
interphalangeal joints. Kydex
plastic forearm and wrist section.
Sizes: SM, MD, LG




#4
Thomas Suspension Splint to dorsiflex the wrist. For
radial palsy. Kydex plastic forearm shell. Outrigger
for thumb. Sizes: SM, MD, LG




#3
Modified Oppenheimer Splint with Coiled Spring Wire
to dorsiflex the wrist for radial palsy. Thumb outrigger.
Hook and loop closure at wrist. Sizes: SM, MD, LG

MORE DETAILED INFORMATION AVAILABLE AT WWW.BUNNELLSPLINTS.COM




#13

Combination of Coiled Spring Wire Oppenheimer to dorsiflex the wrist and Outrigger to extend distal joints (shown) or distal segment (provided) to flex distal joints. Thumb Outrigger. Hook and loop closure at wrist
Sizes: SM, MD, LG



#13A

Combination of Coiled Spring Wire Oppenheimer to dorsiflex the wrist with Reverse Knuckle Bender to extend the metacarpophalangeal joints, and outrigger to extend distal joints. Thumb outrigger. Hook and loop closure at wrist.
Sizes: SM, MD, LG




#13B

Combination of Coiled Spring Wire Oppenheimer to dorsiflex the wrist with Composite Elastic Splint to extend the M.P. and P.I.P. joints. Thumb outrigger. Hook and loop closure at wrist.
Size: SM, MD, LG

MORE DETAILED INFORMATION AVAILABLE AT WWW.BUNNELLSPLINTS.COM

ADJUSTING STATIC POSITIONING

- INCREASE R.O.M.—JOINT STIFFNESS—ELASTIC CONTRACTURES •
- COMFORTABLE PADDING • STAINLESS STEEL CONSTRUCTION •
- EASY TO USE HOOK & LOOP STRAPS TO CONTROL TENSION •



← SIZE	LENGTH OF SPLINT	
	INCHES	CM
XS	3 3/4	9.5
SM	3 7/8	9.8
MD	4 1/8	10.5
LG	4 3/8	11.1
XL	4 1/2	11.4

#18 P.I.P. Extension

Malleable finger platform and M.P. extension for improved function. Provided with optional use straps to secure distal portion of finger.



← SIZE	LENGTH OF SPLINT	
	INCHES	CM
XS	1 1/2	3.6
SM	1 5/8	4
MD	1 3/4	4.3
LG	1 7/8	4.8
XL	2	5.1

#19 P.I.P. and D.I.P. Flexion

Adjustable pads apply pressure to flex distal segment



Positional strap allows for optimal strap placement

#20 M.P. and P.I.P. Flexion with Distal Segment attached

Hook and loop straps to allow positioning and removal of distal segment (above)

M.P. Flexion—Remove distal segment (not shown)

SIZED ACCORDING TO WIDTH ACROSS M.P. JOINT

SIZE	SM	MD	LG
INCHES	2 1/2	3	3 1/2
CM	6.4	7.6	8.9

MORE DETAILED INFORMATION AVAILABLE AT WWW.BUNNELLSPLINTS.COM

LOW PROFILE FINGER SPLINTS



- CLOSED CELL FOAM AND STAINLESS STEEL WIRE CONSTRUCTION
- EASY TO FIT AND ADJUST EXTENSION FORCE
- LATEX FREE
- CONTOURED PADS AND BANDS FOR IMPROVED FIT

SIZING— MEASURE LENGTH FROM DISTAL PALMER CREASE TO CREASE AT D.I.P.
POSITION PROXIMAL BAND OF SPLINT UNDER M.P. HEAD— DISTAL BAND BEHIND D.I.P. JOINT

FW01 PIP EXTENSION

Post injury—joint stiffness


Standard	Length of splint	
FW01-1 XS	2 1/4"	5.7CM
FW01-2 SM	2 5/8"	6.7CM
FW01-3 MD	3"	7.6CM
FW01-4 LG	3 1/2"	8.6CM
FW01-5 XL	4"	10.0CM
Wide	Length of splint	
FW01-1 XS	2 1/4"	5.7CM
FW01-2 SM	2 5/8"	6.7CM
FW01-3 MD	3"	7.6CM
FW01-4 LG	3 1/2"	8.6CM
FW01-5 XL	4"	10.0CM



FW07.....FOR A LITTLE MORE FORCE

Terminal PIP extension

	Length of splint	
FW07-1 XS	2 1/4"	5.7CM
FW07-2 SM	2 5/8"	6.7CM
FW07-3 MD	3"	7.6CM
FW07-4 LG	3 3/8"	8.6CM
FW07-5 XL	3 3/4"	9.5CM




FW02 PIP SPRING EXTENSION—CAPENER TYPE DESIGN

Dynamic extension while allowing active flexion

Joint contractures that are elastic

	Length of splint	
FW02-2 SM	2 1/2"	6.4CM
FW02-3 MD	2 3/4"	7.0CM
FW02-4 LG	3"	7.6CM
FW02-5 XL	3 1/4"	8.3CM



SIZING INFORMATION

BUNNELL SPLINTS FIT RIGHT OR LEFT HAND

EASY TO ADJUST—STAINLESS STEEL CONSTRUCTION

•• *FINGER SPLINTS* ••

SIZED ACCORDING TO LENGTH OF SPLINT

MODEL # 5, 5A			MODEL # 6, 6A Approx. Width at Mid Splint 1"			MODEL # MINI 6, 6A Approx. Width at Mid Splint 5/8"			MODEL # 11, 11A			
SIZE	INCHES	CM	SIZE	INCHES	CM	SIZE	INCHES	CM	SIZE	INCHES	CM	Approx. Width at M.P.
XS	6	15.2	XS	2	5.1	1	1 5/16	3.3	XS	1 1/2	3.6	5/8"
SM	6 1/2	16.5	SM	2 1/2	6.4	2	1 9/16	4	SM	1 5/8	4	7/8"
MD	7	17.8	MD	3	7.6	3	1 7/8	4.8	MD	1 3/4	4.3	1"
LG	7 1/2	19.1	LG	3 1/2	8.9	4	2	5.1	LG	1 7/8	4.8	1 1/8"
XL	8	20.3	XL	4	10.2				XL	2	5.1	1 1/4"

•• *HAND FINGER & WRIST HAND FINGER SPLINTS* ••

SIZED ACCORDING TO WIDTH ACROSS M.P. JOINTS


HAND FINGER SPLINTS

MODEL # 2, 2A, 2B, 9, 9A, 10

WRIST HAND FINGER SPLINTS

MODEL # 1, 1A, 3, 4, 8, 8A, 13, 13A, 13B

SIZE	WIDTH ACROSS M.P. JOINT
XS	2"
SM	2 1/2"
MD	3"
LG	3 1/2"
XL	4"



WENIGER/BUNNEL
TIBURON MEDICAL ENTERPRISES, INC.
915 INDUSTRIAL WAY
SAN JACINTO, CA 92582-3890
PH. (951) 654-2333
EMAIL: SALES@TIBURONMED.COM